

INSTRUCTIONS FOR HOVA-BATOR 1588 INCUBATORS

Cut out instruction panels below and affix to incubator



24 HOURS PRIOR TO SCHEDULED EMBRYOLOGY PRESENTATION

Plug in. Fill the CIRCULATED trough marked 'C' in diagram 1 with warm water.
Calibration of incubator is not needed – no temperature adjustment will be required!




DIAGRAM 1

WHEN EMBRYOLOGY PROJECT IS COMPLETED

Please clean incubator with warm soap and water. Store in original cardboard box.

PLEASE DO NOT LOAN OUT – FOR 3RD GRADE USE ONLY (incubators do not travel well)

Property of: Lancaster County Extension, 444 Cherrycreek Road, Suite A, Lincoln, NE 68528; (402) 441-7180